

AP SERIES

CONTROL AND INTERFACE SYSTEM

For IBERTEST Equipments


To Be Used with:

- Le Chatelier bath
- Automatic mortar mixer
- Jolting table
- Electric flow table
- Abrasion testing machine
- Curing tanks



ADVANTAGES

- **Flexibility:** programmable according to each customer needs and requirements. Allows real time test parameters reading and modification, thanks to a RTC (Real Time Clock) incorporated.
- **Integration:** allowing control and automation implementation.
- **Modularity:** autonomous system, easily replaceable
- **Security:** shortcut or probe disconnection alarm, over tension, component failure, etc.
- **Environment tightness:** designed to be used on industrial environment. IP65 protection level.
- **Robustness:** Minimized maintenance and repair processes. Outputs and inputs isolated to reduce breakdowns and problems identification.
- **Improved interface:** easy and intuitive managing and program.
- **LCD display:** two lines, 16 characters per line. Backlighted screen. Alphanumeric keyboard.
- **Multilanguage:** Spanish, English, French and Portuguese. Other languages can be included on demand.



IBERTEST Internacional S.A.
 European Union Manufacturer
 Ramon y Cajal 35. Daganzo de Arriba,
 Madrid, Spain
 Ph: (+34) 91 884 53 85 Fax: (+34) 91 884 50 02
 sales@ibertestint.com
 www.ibertestint.com